

LDLRAD2 Antibody



PACO01019

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

LDLRAD2

Uniprot

Q5SZI1

Synonyms:

LDLRAD2; Low-density lipoprotein receptor class A domain-containing protein 2

Immunogen:

Synthesized peptide derived from the Internal region of human LDLRAD2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of NIH-3T3 cells using LDLRAD2 Polyclonal Antibody.

