

CACNG7 Antibody



PACO01031

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

CACNG7

Uniprot

P62955

Synonyms:

CACNG7; Voltage-dependent calcium channel gamma-7 subunit; Neuronal voltage-gated calcium channel gamma-7 subunit; Transmembrane AMPAR regulatory protein gamma-7; TARP gamma-7

Immunogen:

Synthesized peptide derived from the C-terminal region of human L-type Ca⁺⁺ CP gamma7.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HepG2 cells using L-type Ca⁺⁺ CP gamma7 Polyclonal Antibody.

