

CHRM4 Antibody



PACO01037

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

CHRM4

Uniprot

P08173

Synonyms:

CHRM4; Muscarinic acetylcholine receptor M4

Immunogen:

Synthesized peptide derived from the Internal region of human mAChR M4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of LOVO cells using mAChR M4 Polyclonal Antibody.

