

# MAST2 Antibody



PACO01045

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300,  
IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

MAST2

**Uniprot**

Q6P0Q8

**Synonyms:**

MAST2; KIAA0807; MAST205; Microtubule-associated serine/threonine-protein kinase 2

**Immunogen:**

Synthesized peptide derived from the Internal region of human MAST205.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of 293 cells using MAST205 Polyclonal Antibody.

