

MPHOSPH9 Antibody



PACO01081

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

MPHOSPH9

Uniprot

Q99550

Synonyms:

MPHOSPH9; MPP9; M-phase phosphoprotein 9


Immunogen:

Synthesized peptide derived from the C-terminal region of human MPP9.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using MPP9 Polyclonal Antibody.