

MRPL24 Antibody



PACO01097

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

MRPL24

Uniprot

Q96A35

Synonyms:

MRPL24; 39S ribosomal protein L24; mitochondrial; L24mt; MRP-L24


Immunogen:

Synthesized peptide derived from the Internal region of human MRP-L24.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc cells using MRP-L24 Polyclonal Antibody.