

MRPL54 Antibody



PACO01104

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

MRPL54

Uniprot

Q6P161

Synonyms:

MRPL54; 39S ribosomal protein L54; mitochondrial; L54mt; MRP-L54


Immunogen:

Synthesized peptide derived from the Internal region of human MRP-L54.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using MRP-L54 Polyclonal Antibody.