

MUTYH Antibody



PACO01123

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

MUTYH

Uniprot

Q9UIF7

Synonyms:

MUTYH; MYH; A/G-specific adenine DNA glycosylase; MutY homolog; hMYH


Immunogen:

Synthesized peptide derived from the Internal region of human MYH.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2 cells using MYH Polyclonal Antibody.