

NOL10 Antibody



PACO01180

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

NOL10

Uniprot

Q9BSC4

Synonyms:

NOL10; Nucleolar protein 10


Immunogen:

Synthesized peptide derived from the C-terminal region of human NOL10.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc cells using NOL10 Polyclonal Antibody.