

NPDC1 Antibody



PACO01185

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

NPDC1

Uniprot

Q9NQX5

Synonyms:

NPDC1; Neural proliferation differentiation and control protein 1; NPDC-1


Immunogen:

Synthesized peptide derived from the Internal region of human NPDC-1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of LOVO cells using NPDC-1 Polyclonal Antibody.