

NRIP3 Antibody



PACO01189

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

NRIP3

Uniprot

Q9NQ35

Synonyms:

NRIP3; C11orf14; Nuclear receptor-interacting protein 3; Sarcoma antigen NY-SAR-105


Immunogen:

Synthesized peptide derived from the C-terminal region of human NRIP3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using NRIP3 Polyclonal Antibody.