

# OR10D4P Antibody



PACO01202

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR10D4P

**Uniprot**

Q8NGN7

**Synonyms:**

OR10D4P; OR10D4; Putative olfactory receptor 10D4

**Immunogen:**


Synthesized peptide derived from the Internal region of human Olfactory receptor 10D4.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of COLO205 cells using Olfactory receptor 10D4 Polyclonal Antibody.