

OR10G6 Antibody



PACO01204

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR10G6

Uniprot

Q8NH81

Synonyms:

OR10G6; OR10G6P; Olfactory receptor 10G6; Olfactory receptor OR11-280

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 10G6.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of COLO205 cells using Olfactory receptor 10G6 Polyclonal Antibody.

