

OR2A7/OR2A4 Antibody



PACO01218

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR2A7/OR2A4

Uniprot

Q96R45/O95047

Synonyms:

OR2A7; Olfactory receptor 2A7; Olfactory receptor OR7-18; OR2A4; OR2A10; Olfactory receptor 2A4; Olfactory receptor 2A10; Olfactory receptor OR6-37

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2A4/7.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HeLa cells using Olfactory receptor 2A4/7 Polyclonal Antibody.

