

OR2J2 Antibody



PACO01221

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR2J2

Uniprot

O76002/Q5SUJ6/Q5SUJ7

Synonyms:

Olfactory receptor 2J2

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2J2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HT29 cells using Olfactory receptor 2J2 Polyclonal Antibody.