

OR2L5 Antibody



PACO01223

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR2L5

Uniprot

Q8NG80

Synonyms:

OR2L5; OR2L11; Olfactory receptor 2L5; Olfactory receptor 2L11; Olfactory receptor OR1-53


Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2L5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using Olfactory receptor 2L5 Polyclonal Antibody.