

OR2M7 Antibody



PACO01224

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR2M7

Uniprot

Q8NG81

Synonyms:

OR2M7; Olfactory receptor 2M7; Olfactory receptor OR1-58

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2M7.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HepG2 cells using Olfactory receptor 2M7 Polyclonal Antibody.

