

OR2T1 Antibody



PACO01225

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR2T1

Uniprot

O43869

Synonyms:

OR2T1; Olfactory receptor 2T1; Olfactory receptor 1-25; OR1-25; Olfactory receptor OR1-61

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2T1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HT29 cells using Olfactory receptor 2T1 Polyclonal Antibody.

