

OR4C12 Antibody



PACO01229

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

OR4C12

Uniprot

Q96R67

Synonyms:

OR4C12; Olfactory receptor 4C12; Olfactory receptor OR11-259


Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 4C12.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using Olfactory receptor 4C12 Polyclonal Antibody.