

# OR4C6 Antibody



PACO01232

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR4C6

**Uniprot**

Q8NH72

**Synonyms:**

OR4C6; Olfactory receptor 4C6; Olfactory receptor OR11-138

**Immunogen:**

Synthesized peptide derived from the Internal region of human Olfactory receptor 4C6.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of HepG2 cells using Olfactory receptor 4C6 Polyclonal Antibody.

