

OR4D1 Antibody



PACO01233

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR4D1

Uniprot

Q15615

Synonyms:

OR4D1; OR4D3; Olfactory receptor 4D1; Olfactory receptor 4D3; Olfactory receptor TPCR16

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 4D1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of Jurkat cells using Olfactory receptor 4D1 Polyclonal Antibody.

