

OR4K14 Antibody



PACO01234

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

OR4K14

Uniprot

Q8NGD5

Synonyms:

OR4K14; Olfactory receptor 4K14; Olfactory receptor OR14-22

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 4K14.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc cells using Olfactory receptor 4K14 Polyclonal Antibody.