

OR52E6 Antibody



PACO01244

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

OR52E6

Uniprot

Q96RD3

Synonyms:

OR52E6; Olfactory receptor 52E6; Olfactory receptor OR11-58

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 52E6.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of MCF7 cells using Olfactory receptor 52E6 Polyclonal Antibody.

