

OR52W1 Antibody



PACO01246

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

OR52W1

Uniprot

Q6IF63

Synonyms:

OR52W1; OR52W1P; Olfactory receptor 52W1; Olfactory receptor OR11-71


Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 52W1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Olfactory receptor 52W1 Polyclonal Antibody.