

# OR56B4 Antibody



PACO01247

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR56B4

**Uniprot**

Q8NH76

**Synonyms:**

OR56B4; Olfactory receptor 56B4; Olfactory receptor OR11-67

**Immunogen:**


Synthesized peptide derived from the C-terminal region of human Olfactory receptor 56B4.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of COLO205 cells using Olfactory receptor 56B4 Polyclonal Antibody.