

OR5AK3P Antibody



PACO01248

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR5AK3P

Uniprot

Q8NH89

Synonyms:

OR5AK3P; OR5AK3; Putative olfactory receptor 5AK3

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 5AK3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of COS7 cells using Olfactory receptor 5AK3 Polyclonal Antibody.

