

# OR5F1 Antibody



PACO01250

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR5F1

**Uniprot**

O95221

**Synonyms:**

OR5F1; Olfactory receptor 5F1; Olfactory receptor 11-10; OR11-10; Olfactory receptor OR11-167

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 5F1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of Jurkat cells using Olfactory receptor 5F1 Polyclonal Antibody.

