

OR5K1 Antibody



PACO01251

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR5K1

Uniprot

Q8NHB7

Synonyms:

OR5K1; Olfactory receptor 5K1; HTPCRX10; Olfactory receptor OR3-8


Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 5K1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using Olfactory receptor 5K1 Polyclonal Antibody.