

OR5P3 Antibody



PACO01253

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

OR5P3

Uniprot

Q8WZ94

Synonyms:

OR5P3; Olfactory receptor 5P3; Olfactory receptor OR11-94; Olfactory receptor-like protein JCG1


Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 5P3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using Olfactory receptor 5P3 Polyclonal Antibody.