

OR6C3 Antibody



PACO01254

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

OR6C3

Uniprot

Q9NZP0

Synonyms:

OR6C3; Olfactory receptor 6C3; HSA8

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 6C3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HepG2 cells using Olfactory receptor 6C3 Polyclonal Antibody.

