

# OR6C68 Antibody



PACO01255

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

OR6C68

**Uniprot**

A6NDL8

**Synonyms:**

OR6C68; Olfactory receptor 6C68

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 6C68.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of HuvEc cells using Olfactory receptor 6C68 Polyclonal Antibody.

