

OR6K2 Antibody



PACO01256

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

OR6K2

Uniprot

Q8NGY2

Synonyms:

OR6K2; Olfactory receptor 6K2; Olfactory receptor OR1-17


Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 6K2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Olfactory receptor 6K2 Polyclonal Antibody.