

OR6T1 Antibody



PACO01257

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

OR6T1

Uniprot

Q8NGN1

Synonyms:

OR6T1; Olfactory receptor 6T1; Olfactory receptor OR11-277


Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 6T1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Olfactory receptor 6T1 Polyclonal Antibody.