

OR7C2 Antibody



PACO01260

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

OR7C2

Uniprot

O60412

Synonyms:

OR7C2; OR7C3; Olfactory receptor 7C2; Olfactory receptor 19-18; OR19-18; Olfactory receptor 7C3; Olfactory receptor OR19-22

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 7C2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HuvEc cells using Olfactory receptor 7C2 Polyclonal Antibody.

