

OR8S1 Antibody



PACO01263

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

OR8S1

Uniprot

Q8NH09

Synonyms:

OR8S1; Olfactory receptor 8S1


Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 8S1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HT29 cells using Olfactory receptor 8S1 Polyclonal Antibody.