

CDKN2B Antibody



PACO01277

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

CDKN2B

Uniprot

P42772

Synonyms:

CDKN2B; MTS2; Cyclin-dependent kinase 4 inhibitor B; Multiple tumor suppressor 2;
MTS-2; p14-INK4b; p15-INK4b; p15INK4B


Immunogen:

Synthesized peptide derived from the C-terminal region of human p15.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using p15 Polyclonal Antibody.