

PDRG1 Antibody



PACO01312

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

PDRG1

Uniprot

Q9NUG6

Synonyms:

PDRG1; C20orf126; PDRG; p53 and DNA damage-regulated protein 1


Immunogen:

Synthesized peptide derived from the Internal region of human PDRG1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2 cells using PDRG1 Polyclonal Antibody.