

PDXDC1 Antibody



PACO01313

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

PDXDC1

Uniprot

Q6P996

Synonyms:

PDXDC1; KIAA0251; Pyridoxal-dependent decarboxylase domain-containing protein 1


Immunogen:

Synthesized peptide derived from the C-terminal region of human PDXDC1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc cells using PDXDC1 Polyclonal Antibody.