

PHF3 Antibody



PACO01330

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

PHF3

Uniprot

Q92576

Synonyms:

PHF3; KIAA0244; PHD finger protein 3


Immunogen:

Synthesized peptide derived from the C-terminal region of human PHF3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using PHF3 Polyclonal Antibody.