

# PPP1R14D Antibody



PACO01358

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

PPP1R14D

**Uniprot**

Q9NXH3

**Synonyms:**

PPP1R14D; GBPI; Protein phosphatase 1 regulatory subunit 14D; Gastrointestinal and brain-specific PP1-inhibitory protein 1; GBPI-1

**Immunogen:**

Synthesized peptide derived from the Internal region of human PPP1R14D.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of COLO205 cells using PPP1R14D Polyclonal Antibody.

