

# PSKH1 Antibody



PACO01370

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

PSKH1

**Uniprot**

P11801

**Synonyms:**

PSKH1; Serine/threonine-protein kinase H1; Protein serine kinase H1; PSK-H1

**Immunogen:**

Synthesized peptide derived from the Internal region of human PSK-H1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of A549 cells using PSK-H1 Polyclonal Antibody.

