

# RAB35 Antibody



PACO01384

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:****Gene ID:**

RAB35

**Uniprot**

Q15286

**Synonyms:**

RAB35; RAB1C; RAY; Ras-related protein Rab-35; GTP-binding protein RAY; Ras-related protein Rab-1C

**Immunogen:**

Synthesized peptide derived from the Internal region of human Rab 35.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of 293 cells using Rab 35 Polyclonal Antibody.

