

# RABEP2 Antibody



PACO01388

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

RABEP2

**Uniprot**

Q9H5N1

**Synonyms:**

RABEP2; RABPT5B; Rab GTPase-binding effector protein 2; Rabaptin-5beta

**Immunogen:**

Synthesized peptide derived from the Internal region of human Rabaptin-5Beta.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of HuvEc cells using Rabaptin-5Beta Polyclonal Antibody.

