

RPL18 Antibody



PACO01437

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

RPL18

Uniprot

Q07020

Synonyms:

RPL18; 60S ribosomal protein L18


Immunogen:

Synthesized peptide derived from the Internal region of human Ribosomal Protein L18.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using Ribosomal Protein L18 Polyclonal Antibody.