

RPL27A Antibody



PACO01439

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

RPL27A

Uniprot

P46776

Synonyms:

RPL27A; 60S ribosomal protein L27a


Immunogen:

Synthesized peptide derived from the C-terminal region of human Ribosomal Protein L27A.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using Ribosomal Protein L27A Polyclonal Antibody.