

RPL30 Antibody



PACO01442

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

RPL30

Uniprot

P62888

Synonyms:

RPL30; 60S ribosomal protein L30


Immunogen:

Synthesized peptide derived from the C-terminal region of human Ribosomal Protein L30.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAW264.7 cells using Ribosomal Protein L30 Polyclonal Antibody.