

SCN4B Antibody



PACO01477

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

SCN4B

Uniprot

Q8IWT1

Synonyms:

SCN4B; Sodium channel subunit beta-4


Immunogen:

Synthesized peptide derived from the Internal region of human Scn4b.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using Scn4b Polyclonal Antibody.