

SUPT3H Antibody



PACO01532

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

SUPT3H

Uniprot

O75486

Synonyms:

SUPT3H; SPT3; Transcription initiation protein SPT3 homolog; SPT3-like protein

Immunogen:

Synthesized peptide derived from the Internal region of human SPT3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HepG2 cells using SPT3 Polyclonal Antibody.

