

# HTR1E Antibody



PACO01534

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

HTR1E

**Uniprot**

P28566

**Synonyms:**

HTR1E; 5-hydroxytryptamine receptor 1E; 5-HT-1E; 5-HT1E; S31; Serotonin receptor 1E

**Immunogen:**

Synthesized peptide derived from the Internal region of human SR-1E.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of HeLa cells using SR-1E Polyclonal Antibody.

