

# STAM2 Antibody



PACO01546

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300,  
IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

STAM2

**Uniprot**

O75886

**Synonyms:**

STAM2; HBP; Signal transducing adapter molecule 2; STAM-2; Hrs-binding protein

**Immunogen:**

Synthesized peptide derived from human STAM2 around the non-phosphorylation site of Y192.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of NIH-3T3 cells using STAM2 Polyclonal Antibody.

